




# XUV Phase Retarder AURORA

Our XUV Phase Retarder acts as a quarter waveplate to turn linearly polarized XUV light into circularly polarized light without introducing noticeable dispersion. The phase retarder achieves close-to-circular polarization of  $P_c=0.75$  and obtains  $> 25\%$  transmission around 66 eV photon energy, where the Ni  $M_2/M_3$  edge locates [1]. A broad spectral range from 40 to 85 eV is supported to cover the  $M_2/M_3$  edge of 3 transition metals - Fe, Co and Ni.

The retarder uses a transmission optimized, four mirror grazing incidence reflection geometry that induces a quarter wave phase offset between the s- and p-polarization components of a linearly polarized input XUV beam. At the Ni  $M_2/M_3$  transition up to 3% dichroism contrast is observed (i.e.  $> 85\%$  of the theoretical value) [1].

A clear aperture of 3 mm will allow the low divergent XUV to pass through without clipping. Our XUV

Phase Retarder is ideally geared to be combined with ultrafast high-harmonic XUV sources adding spin-sensitivity to conventional laser based pump-probe experiments via attosecond magnetic circular dichroism detection.



### Key Product Features:


- Create circularly polarized XUV light without adding dispersion (ideal for attosecond applications)
- Up to  $\approx 40\%$  max transmission
- Broad spectral range covering 40 – 85 eV
- Compact design, with minimum change in the delay line
- Entrance and exit iris diaphragms for alignment
- Easy to implement in existing setups

### Characteristics:


Spectral Range	40-85 eV
Avg Transmission	$>25\%$
Max Ellipticity $P_c^1$	0.75 @ 66 eV
Extra Beam Path	3 mm
Clear Aperture	3 mm
Footprint	40mm (D) x 92mm (L)

<sup>1</sup>Degree of Ellipticity expressed as Stokes Parameters  $P_c=S3/S0$





Phase and Transmission, p-Pol



Phase and Transmission, s-Pol

Element	Fe	Co		Ni	
Edge	$M_2/M_3$	$M_2$	$M_3$	$M_2$	$M_3$
Energy	52.7 eV	58.9 eV	59.9 eV	66.2 eV	68 eV
Transmission	9%	17%	18%	26%	30%

### References:

- [1] Siegrist, F., Gessner, J.A., Ossiander, M. et al. Light-wave dynamic control of magnetism. *Nature* **571**, 240-244(2019).