



**UltraFast**  
Innovations


YOUR KEY to innovation and success



## Cavity-Ringdown (CRD) Reflectometer and Loss Meter

### GLACIER®

Our reflectometer GLACIER uses the extreme sensitivity of cavity ring-down spectroscopy to quantify the losses of advanced optical coatings down to 5 ppm. As a typical application the device can characterize supra-mirrors with up to 99.9995% reflectivity. Conventional absorption and reflection measurements are not sufficiently sensitive to quantify today's super-reflective mirror coatings and are typically limited to the  $>1000$  ppm range (corresponding to  $<99.9\%$  reflectivity). Cavity ring-down spectroscopy measures optical losses by the decay of the energy stored inside a cavity.



The technique reaches unrivalled sensitivity because the losses are experienced over and over again after every round trip inside the cavity. Smaller losses lead to longer intra-cavity dwell time thereby automatically increasing the mea-

surement precision. The device features high-speed data acquisition and allows to record measurements within seconds. It is delivered complete with a computer and a user-friendly software interface for acquisition and real-time analysis.

### Key Product Features:

- Reflectivity measurements
  - Reflectivities up to 99.9995%
  - Various angles of incidence:  $5^{\circ}$ - $45^{\circ}$  (and  $0^{\circ}$ )
  - s and p polarization (separately)
- Antireflective coating characterization  
Reflectivity down to 0.0005% (5 ppm)
- Simple and reproducible alignment for 0.5", 1" and 2" optics
- Spring-loaded mirror fixtures for reproducible mounting without strain
- Computer and user-friendly software interface included
- High-speed data acquisition and real-time analysis  
Complete measurement and analysis within seconds
- Available wavelengths: 375-1550 nm  
other wavelengths upon request
- Standard Footprint: 90 x 45 cm<sup>2</sup>

UltraFast Innovations GmbH  
Dieselstr. 5  
85748 Garching  
Germany


phone: +49 89 36039 - 437  
fax: +49 89 36039 - 453  
info@ultrafast-innovations.com  
www.ultrafast-innovations.com



	GLACIER	GLACIER*	GLACIER**
Number of wavelengths	One	Two	Three
Wavelength range	375-1550 nm	375-1550 nm	375-1550 nm
Footprint	90 x 45 cm <sup>2</sup>	90 x 55 cm <sup>2</sup>	

### Working Principle:

GLACIER uses the principle of reflectivity/loss measurements with cavity ring-down spectroscopy based on very low losses at each mirror bounce. The laser pulses travel inside a cavity experiencing optical losses over and over again during each round trip.



Sketch of the working principle of GLACIER, measurement and fitting procedure.

The device measures the time-dependent intensity  $I(t)$  leaked through an end mirror of the cavity (center). The signal decays with a time constant depending on the intra-cavity losses and can be fitted to the following exponential function:


$$I(t) = I(t_0) \cdot \exp\left(-\frac{t}{\tau}\right)$$

The time constant  $\tau$  is inversely proportional to the optical losses  $(1-R)$  of the cavity with total reflectivity  $R$ :

$$\tau = \frac{n}{c} \cdot \frac{l}{(1-R)}$$

where  $n$  is the refractive index,  $c$  is the speed of light, and  $l$  is the cavity length.

### Sample Measurement:



Typical GLACIER measurement of an ultra-high-reflective mirror for 1030 nm. To obtain the data, the cavity losses with and without the sample were measured and subtracted. This provides an absolute measurement of the test mirror.